ENGLISH

Product Description

REGENN I.V.[®] hydrocolloid dressings are sterile dressings that are designed to help protect peripheral intravenous and central catheter insertion sites and assist in securing catheters in place.

The REGENN I.V.[®] dressings have a circular construction that consist of a skin contact surface layer of thin hydrocolloid adhesive with an embossed outer ring of a multi-layer hydrocolloid adhesive that is design to provide secure adhesion at the periphery of the dressing. Along the outer perimeter of the dressing, there is a notch in the multi layered hydrocolloid adhesive that is intended to accommodate the egress of a peripheral or central line catheter.

The top layer of the dressing consists of a low friction, breathable, polyurethane film.

A "Z-fold" designed silicone-release paper system allows for ease of application and minimizes the risk of gloves sticking to the adhesive surface of the dressing. REGENN I.V.[®] dressing's thin hydrocolloid adhesive is gentle even on sensitive skin, can absorb minimal exudate, and allows the dressing to adhere to dry and moist skin surface. The transparency of the dressing allows for visualization of the insertion site without the need for removal.

The barrier properties of the dressing protect against external contaminants and allow patients to shower and the dressing's enhanced flexibility and conformability help the dressing to adhere well in difficult anatomical locations.

As long as dressing remains secure to the skin, the dressing provides a viral barrier.

Product Indications: REGENN I.V.[®] dressings can be used to cover and protect catheter sites and to secure devices to the skin REGENN I.V.[®] dressing's common applications include:

- Peripheral I.V. catheters.
- Other intravascular catheters and percutaneous devices.

Product Contraindications: Known sensitivity to the dressing or any of its components.

Precautions

- 1 - Stop any bleeding at the site before applying the dressing.
- 2 - Do not stretch the dressing during application as this could impact adhesive performance.
- 3 - Make sure the skin is dry and free of soap residue and lotion to ensure good adhesion.
- 4 - If infection should occur that is related to the insertion site, appropriate medical treatment should be initiated. This dressing may continue to be used under appropriate medical supervision.

Site preparation

- 1 - Prepare the site according to institutional protocol.
- 2 - Clipping of hair at the site may improve dressing adherence. Shaving is not recommended.
- 3 - Preps and protectants are normally not needed

Product Application

- 1 - Cleanse and dry around insertion site as indicated by standard protocol.
- 2 - Remove dressing from packaging.
- 3 - With one hand, grasp white release paper at end with arrows and peel back part way.
- 4 - With other hand, grasp the clear film at end with arrows and turn dressing over so that notch in dressing faces direction of catheter egress from insertion site.
- 5 - Press dressing to the skin over the insertion site.
- 6 - Remove white release paper and smooth dressing to skin.
- 7 - Remove clear film from other side of dressing and smooth dressing to the skin.
- 8 - Pinch dressing around the catheter at exit notch to help secure tube to skin.
- 9 - Apply tape strip for additional tube support if needed.
- 10 - Writing on the dressing is not recommended.

Site Care

- 1 - The site should be observed daily for signs of infection or other complications. If infection is suspected, remove the dressing, inspect the site directly, and determine appropriate medical intervention. Infection may be signaled by fever, pain, redness swelling, or unusual odor or discharge.
- 2 - Change the dressing according to institution protocol, or when the barrier properties have been compromised.
- 3 - Dressing can remain in place for up to 7 days.
- 4 - If leakage of exudates should occur the dressing should be changed immediately.

Product Removal

- 1 - With one hand, lift edge of dressing gently from the skin.
- 2 - Use the other hand to press down on skin to facilitate removal.
- 3 - Work your way around the dressing using this technique.
- 4 - Roll dressing off skin without pulling on catheter tubing or the catheter hub.
- 5 - Dispose of used dressing as per local regulation.

Shelf Life and Storage Conditions

- Store in a cool dry place.
- Do not use if primary package is visibly damage or broken.
- Do not Re-sterilize or reuse, method of sterilization by radiation.
- Three (3) year shelf life
- Latex free

Laboratory testing has demonstration that REGENN I.V.[®] dressings provide a viral barrier (HIV-1 and HBV) while the dressing remains intact without separation from the skin. The use of this dressing neither guarantees nor warrants protection against AIDS or Hepatitis B Virus transmission).

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